
















Set Code:	01.037.009E
Set Barcode:	 (01)10886982273543
Description:	TFNA 11 X 130 DEG NAIL EV
Case Report ID:	

Surgery Date and Time:	
Surgeon:	
Operating Room:	
Case ID:	

Primary Tray

	Product #	Product Desc.	QTY	UDI-DI	LOT Numbers	UDI Barcode
1	NA0525	TFNA SOFT CASE FOR STERILE NAILS BLADES AND SCREWS	1	10886982260604		 (01)10886982260604
2	04.037.142S	11MM/130 DEG TI CANN TFNA 170MM - STERILE	1	10886982096548		 (01)10886982096548
3	04.037.143S	11MM/130 DEG TI CANN TFNA 200MM - STERILE	1	10886982096555		 (01)10886982096555
4	04.037.144S	11MM/130 DEG TI CANN TFNA 235MM/RIGHT - STERILE	1	10886982096562		 (01)10886982096562
5	04.037.145S	11MM/130 DEG TI CANN TFNA 235MM/LEFT - STERILE	1	10886982096579		 (01)10886982096579
6	04.037.150S	11MM/130 DEG TI CANN TFNA 300MM/RIGHT - STERILE	1	10886982096586		 (01)10886982096586
7	04.037.151S	11MM/130 DEG TI CANN TFNA 300MM/LEFT - STERILE	1	10886982096593		 (01)10886982096593
8	04.037.152S	11MM/130 DEG TI CANN TFNA 320MM/RIGHT - STERILE	1	10886982096609		 (01)10886982096609
9	04.037.153S	11MM/130 DEG TI CANN TFNA 320MM/LEFT - STERILE	1	10886982096616		 (01)10886982096616
10	04.037.154S	11MM/130 DEG TI CANN TFNA 340MM/RIGHT - STERILE	1	10886982096623		 (01)10886982096623
11	04.037.155S	11MM/130 DEG TI CANN TFNA 340MM/LEFT - STERILE	1	10886982096630		 (01)10886982096630
12	04.037.156S	11MM/130 DEG TI CANN TFNA 360MM/RIGHT - STERILE	1	10886982096647		 (01)10886982096647
13	04.037.157S	11MM/130 DEG TI CANN TFNA 360MM/LEFT - STERILE	1	10886982096654		 (01)10886982096654

Set Code:	01.037.009E
Set Barcode:	 (01)10886982273543
Description:	TFNA 11 X 130 DEG NAIL EV
Case Report ID:	

Surgery Date and Time:	
Surgeon:	
Operating Room:	
Case ID:	

Primary Tray

	Product #	Product Desc.	QTY	UDI-DI	LOT Numbers	UDI Barcode
14	04.037.158S	11MM/130 DEG TI CANN TFNA 380MM/RIGHT - STERILE	1	10886982096661		 (01)10886982096661
15	04.037.159S	11MM/130 DEG TI CANN TFNA 380MM/LEFT - STERILE	1	10886982096678		 (01)10886982096678
16	04.037.160S	11MM/130 DEG TI CANN TFNA 400MM/RIGHT - STERILE	1	10886982096685		 (01)10886982096685
17	04.037.161S	11MM/130 DEG TI CANN TFNA 400MM/LEFT - STERILE	1	10886982096692		 (01)10886982096692
18	04.037.162S	11MM/130 DEG TI CANN TFNA 420MM/RIGHT - STERILE	1	10886982096708		 (01)10886982096708
19	04.037.163S	11MM/130 DEG TI CANN TFNA 420MM/LEFT - STERILE	1	10886982096715		 (01)10886982096715
20	04.037.164S	11MM/130 DEG TI CANN TFNA 440MM/RIGHT - STERILE	1	10886982096722		 (01)10886982096722
21	04.037.165S	11MM/130 DEG TI CANN TFNA 440MM/LEFT - STERILE	1	10886982096739		 (01)10886982096739
22	04.037.166S	11MM/130 DEG TI CANN TFNA 460MM/RIGHT - STERILE	1	10886982096746		 (01)10886982096746
23	04.037.167S	11MM/130 DEG TI CANN TFNA 460MM/LEFT - STERILE	1	10886982096753		 (01)10886982096753
24	04.037.168S	11MM/130 DEG TI CANN TFNA 480MM/RIGHT - STERILE	1	10886982096760		 (01)10886982096760
25	04.037.169S	11MM/130 DEG TI CANN TFNA 480MM/LEFT - STERILE	1	10886982096777		 (01)10886982096777