
















Set Code:	01.037.014E
Set Barcode:	 (01)10886982273581
Description:	TFNA 11 X 125 DEG NAIL EV
Case Report ID:	

Surgery Date and Time:	
Surgeon:	
Operating Room:	
Case ID:	

**Primary Tray**

	Product #	Product Desc.	QTY	UDI-DI	LOT Numbers	UDI Barcode
1	NA0525	TFNA SOFT CASE FOR STERILE NAILS BLADES AND SCREWS	1	10886982260604		 (01)10886982260604
2	04.037.112S	11MM/125 DEG TI CANN TFNA 170MM - STERILE	1	10886982096302		 (01)10886982096302
3	04.037.113S	11MM/125 DEG TI CANN TFNA 200MM - STERILE	1	10886982096319		 (01)10886982096319
4	04.037.114S	11MM/125 DEG TI CANN TFNA 235MM/RIGHT - STERILE	1	10886982096326		 (01)10886982096326
5	04.037.115S	11MM/125 DEG TI CANN TFNA 235MM/LEFT - STERILE	1	10886982096333		 (01)10886982096333
6	04.037.120S	11MM/125 DEG TI CANN TFNA 300MM/RIGHT - STERILE	1	10886982096340		 (01)10886982096340
7	04.037.121S	11MM/125 DEG TI CANN TFNA 300MM/LEFT - STERILE	1	10886982096357		 (01)10886982096357
8	04.037.122S	11MM/125 DEG TI CANN TFNA 320MM/RIGHT - STERILE	1	10886982096364		 (01)10886982096364
9	04.037.123S	11MM/125 DEG TI CANN TFNA 320MM/LEFT - STERILE	1	10886982096371		 (01)10886982096371
10	04.037.124S	11MM/125 DEG TI CANN TFNA 340MM/RIGHT - STERILE	1	10886982096388		 (01)10886982096388
11	04.037.125S	11MM/125 DEG TI CANN TFNA 340MM/LEFT - STERILE	1	10886982096395		 (01)10886982096395
12	04.037.126S	11MM/125 DEG TI CANN TFNA 360MM/RIGHT - STERILE	1	10886982096401		 (01)10886982096401
13	04.037.127S	11MM/125 DEG TI CANN TFNA 360MM/LEFT - STERILE	1	10886982096418		 (01)10886982096418

Set Code:	01.037.014E
Set Barcode:	 (01)10886982273581
Description:	TFNA 11 X 125 DEG NAIL EV
Case Report ID:	

Surgery Date and Time:	
Surgeon:	
Operating Room:	
Case ID:	

**Primary Tray**

	Product #	Product Desc.	QTY	UDI-DI	LOT Numbers	UDI Barcode
14	04.037.128S	11MM/125 DEG TI CANN TFNA 380MM/RIGHT - STERILE	1	10886982096425		 (01)10886982096425
15	04.037.129S	11MM/125 DEG TI CANN TFNA 380MM/LEFT - STERILE	1	10886982096432		 (01)10886982096432
16	04.037.130S	11MM/125 DEG TI CANN TFNA 400MM/RIGHT - STERILE	1	10886982096449		 (01)10886982096449
17	04.037.131S	11MM/125 DEG TI CANN TFNA 400MM/LEFT - STERILE	1	10886982096456		 (01)10886982096456
18	04.037.132S	11MM/125 DEG TI CANN TFNA 420MM/RIGHT - STERILE	1	10886982096463		 (01)10886982096463
19	04.037.133S	11MM/125 DEG TI CANN TFNA 420MM/LEFT - STERILE	1	10886982096470		 (01)10886982096470
20	04.037.134S	11MM/125 DEG TI CANN TFNA 440MM/RIGHT - STERILE	1	10886982096487		 (01)10886982096487
21	04.037.135S	11MM/125 DEG TI CANN TFNA 440MM/LEFT - STERILE	1	10886982096494		 (01)10886982096494
22	04.037.136S	11MM/125 DEG TI CANN TFNA 460MM/RIGHT - STERILE	1	10886982096500		 (01)10886982096500
23	04.037.137S	11MM/125 DEG TI CANN TFNA 460MM/LEFT - STERILE	1	10886982096517		 (01)10886982096517
24	04.037.138S	11MM/125 DEG TI CANN TFNA 480MM/RIGHT - STERILE	1	10886982096524		 (01)10886982096524
25	04.037.139S	11MM/125 DEG TI CANN TFNA 480MM/LEFT - STERILE	1	10886982096531		 (01)10886982096531